

Sectional Phantom Series (Plain X-rays)

Sectional Phantoms allow for imaging of individual anatomy as needed.



Head



41926-000 Head (Opaque)



41926-010 Head (Transparent)

Optional parts:
Adjustable head supporter

Stand-alone design can be used with the adjustable head positioning stand to demonstrate accurate skull positioning.

* Opaque model shown.



Body



41926-060 Thorax (Opaque)

41926-070 Thorax (Transparent)

Includes thoracic skeletal system with embedded mediastinal space and bronchus to provide realistic imaging. The scapulae are rotated outside of the lung fields for proper PA chest imaging.

* Opaque model shown.



41926-080 Pelvic (Opaque)

Includes lumbar/sacral spine, pelvic bony anatomy and proximal femurs.





Arm




41926-140
Right Elbow
(Opaque)

Normal flexion range allows for AP/lateral and partial flexion views with one phantom.

* Transparent model shown.

41926-150
Right Elbow
(Transparent)



41926-020
Right Hand
(Opaque)

41926-030
Right Hand
(Transparent)

* Transparent model shown.



41926-040
Left Hand
(Opaque)

41926-050
Left Hand
(Transparent)

* Opaque model shown.



Leg




41926-180
Right Knee
(Opaque)

Freely movable patella and joint allows for realistic positioning of the knee for AP/lateral, oblique, sunrise and tunnel views.

* Transparent model shown.

41926-190
Right Knee
(Transparent)



41926-100
Right Foot
(Opaque)

41926-110
Right Foot
(Transparent)

* Transparent model shown.



41926-120
Left Foot
(Opaque)

41926-130
Left Foot
(Transparent)

* Opaque model shown.

