

# Featherweight® Model FWSTSII

## Single Short Endoscopy Cart

*Innovation Designed Around You*

**PHS  
WEST  
INC.**

### Technical Specifications

When the requirement is for something small, simple, and cost effective, our Featherweight® Single Short Endoscopy Cart fits the need.



### Ergo-Express® Model FWSTSII | Single Short Endoscopy Carts

#### *Specifications*

- Model Number: FWSTSII
- Dimensions:
  - Depth: 24"
  - Width: 25"
  - Height: 5=45.5" | 8=48.5"
  - Weight: Dependent upon cart size
  - Shipping Weight: Dependent upon cart size
- Extrusion Dimensions and Materials
  - 30mm x 60mm Aluminum
- Caster Wheels
  - Plated Carbon Steel
- Bumpers
  - Vinyl
- Drawer Options
  - Sizes: 4" | 6" | 8" | 12"
  - Width: 20.25"
  - Height: 3.5" | 5.25" | 7.125" | 10.75"
  - Depth: 18.375"
- Materials
  - Plastic Panels
    - ABS
  - Shelves and Deck
    - Aluminum
  - Work Surface
    - PVC
- Three Year Limited Warranty
- [24/7-365 Technical and Operational Support for the life of the product](#)

Additional technical and operational specifications are available for custom designs.

