

# Featherweight® Model 2000RFWSTSII Motorized Single Short Endoscopy Cart

Innovation Designed Around You

**PHS  
WEST  
INC.**

## Technical Specifications

When the requirement is for something small, simple, and cost effective, our Featherweight® Motorized Single Short Endoscopy Cart fits the need.



## Ergo-Express® Model 2000RFWSTSII | Motorized Single Short Endoscopy Carts

### Specifications

- Model Number: 2000RFWSTSII
- Dimensions:
  - Depth: 24"
  - Width: 32.5"
  - Height: 48.5"
  - Weight: Dependent upon cart size
  - Shipping Weight: Dependent upon cart size
- Extrusion Dimensions and Materials
  - 30mm x 60mm Aluminum
- Caster Wheels
  - Plated Carbon Steel
- Bumpers
  - Vinyl
- Drawer Options
  - Sizes: 4" | 6" | 8" | 12"
  - Width: 20.25"
  - Height: 3.5" | 5.25" | 7.125" | 10.75"
  - Depth: 18.375"
- Materials
  - Plastic Panels
    - ABS
  - Shelves and Deck
    - Aluminum
  - Work Surface
    - PVC
- ETL/ETLC Certified: UL 583:2012
  - Ed. 10 +R:25May2016, ULC C583-1974/01/01
- Controller: Microprocessor-based
- Charging System
  - 120VAC
  - 2.4 A
  - Watts: 288
  - Electrical Phase: 1
  - Charging Cord: (UL) 384516 SJT 3X14AWG 105°C 300V VW-1
  - Plug Type: 5-15P
  - Drive System Voltage: 24VDC
- Battery Specifications
  - 35 AH or 55 AH Option
  - AGM Battery: Stand-By, UPS, High-Rate, Telecomm
  - Absorbed Glass Mat (AGM) - Battery
  - Synonyms: SLA, VRLA, AGM, Absorbed Glass Mat
  - Sealed Recombinant
  - Dot Description: Wet Battery, Non-Spillable
  - Chemical Name: Sealed Lead Acid, Secondary Battery
- Drive System
  - 24 VDC, 2500 RPM, 5.00 A Maximum, PMDC 4 Brush Motor
  - Electromechanical Holding Brake, rated at 50lb./in. (5.7 Nm), with override and controller feedback
  - Top Speed: 3.00 MPH
- Three Year Limited Warranty
- [24/7-365 Technical and Operational Support for the life of the product](#)

Additional technical and operational specifications are available for custom designs.

