

## IV Prep Cart (SP21W4, SP24W5, SP27W6, SP30W6, and SP33W7 with SPIVP)


Keep pace with supply storage and transport. Keep pace with work flow. Keep pace with patient care.

Pace carts are mobile, easy to maneuver, secure, and adaptable to change. The IV Prep Cart includes an IV pole, upper accessory rail, set of 4 tilt bins, locking Sharps container, and glove box. A variety of drawer configurations and exterior and drawer accessories lets you target supply storage for specific treatments and departments.



Product Shown: IV Prep Cart (SP27W6 with SPIVP)

## IV Prep Cart (SP21W4, SP24W5, SP27W6, SP30W6, and SP33W7 with SPIVP)



### Materials

Shell: **HDPE core with aluminum exterior**  
Base: **16-gauge steel with powder-coat finish**  
Bumpers: **ABS thermoplastic polymer**  
Corner extrusions: **Anodized aluminum**  
Work surface: **Antimicrobial ABS thermoplastic polymer**  
Drawers: **18-gauge steel with powder-coat finish**  
Drawer pulls: **Antimicrobial ABS thermoplastic polymer**  
Drawer fronts: **ABS thermoplastic polymer**  
Drawer glides: **Zinc plated, steel ball bearing**  
Casters: **Heavy duty, 5" full-swivel, ball bearing with steel plate and frame**

### Dimensions (nominal)

#### Pace 21

Exterior: **29.25"w x 21.25"d x 33"h**  
Weight: **110 lbs**

#### Pace 24

Exterior: **29.25"w x 21.25"d x 36"h**  
Weight: **120 lbs**

#### Pace 27

Exterior: **29.25"w x 21.25"d x 39.25"h**  
Weight: **135 lbs**

#### Pace 30

Exterior: **29.25"w x 21.25"d x 42.5"h**  
Weight: **145 lbs**

#### Pace 33

Exterior: **29.25"w x 21.25"d x 45.5"h**  
Weight: **155 lbs**

### Options

- 5 case colors and 11 drawer colors
- Keyed, keypad, proximity, and dual access locks
- Back- and side-mounted accessories
- Drawer organization accessories
- IV Prep Cart configuration (SPIVP) can be ordered for any standard size Pace cart



*Storage solutions for healthcare*

1239 Comstock Street | Marne, Michigan 49435  
888-435-2256 phone | 616-743-5944 fax | [innerspacehealthcare.com](http://innerspacehealthcare.com)