

## Pace Pharmacy Cart (SP21W4, SP24W5, SP27W6, SP30W6, and SP33W7 with SPPVP)


Keep pace with supply storage and transport. Keep pace with work flow. Keep pace with patient care.

Pace carts are mobile, easy to maneuver, secure, and adaptable to change. The Pharmacy Cart has a glove box holder, flip-up work surface, chart holder, disinfectant holder, push handle, 2 security lock boxes, and 2 3"-high pharmacy trays. A variety of drawer configurations and exterior and drawer accessories lets you target supply storage for specific treatments and departments.



Product Shown: Pharmacy Cart (SP27W6 with SPPVP)

# Pace Pharmacy Cart (SP21W4, SP24W5, SP27W6, SP30W6, and SP33W7 with SPPVP)



## Materials

Shell: **HDPE core with aluminum exterior**  
Base: **16-gauge steel with powder-coat finish**  
Bumpers: **ABS thermoplastic polymer**  
Corner extrusions: **Anodized aluminum**  
Work surface: **Antimicrobial ABS thermoplastic polymer**  
Drawers: **18-gauge steel with powder-coat finish**  
Drawer pulls: **Antimicrobial ABS thermoplastic polymer**  
Drawer fronts: **ABS thermoplastic polymer**  
Drawer glides: **Zinc plated, steel ball bearing**  
Casters: **Heavy duty, 5" full-swivel, ball bearing with steel plate and frame**

## Dimensions (nominal)

### Pace 21

Exterior: **29.25"w x 21.25"d x 33"h**  
Weight: **110 lbs**

### Pace 24

Exterior: **29.25"w x 21.25"d x 36"h**  
Weight: **120 lbs**

### Pace 27

Exterior: **29.25"w x 21.25"d x 39.25"h**  
Weight: **135 lbs**

### Pace 30

Exterior: **29.25"w x 21.25"d x 42.5"h**  
Weight: **145 lbs**

### Pace 33

Exterior: **29.25"w x 21.25"d x 45.5"h**  
Weight: **155 lbs**

## Options

- 5 case colors and 11 drawer colors
- Keyed, keypad, proximity, and dual access locks
- Back- and side-mounted accessories
- Drawer organization accessories
- Pharmacy Cart configuration (SPPVP) can be ordered for any standard size Pace cart