

Pace Transport Cart (SP21WC, SP24WC, SP27WC, SP30WC, and SP33WC)


Keep pace with supply storage and transport. Keep pace with work flow. Keep pace with patient care.

Pace carts are mobile, easy to maneuver, secure, and adaptable to change. The FlexCell Cart has an open front and full FlexCell interior. It has a white case. A variety of exterior accessories lets you target supply storage for specific treatments and departments.



Product Shown: FlexCell Cart (SP27WC)

Pace Transport Cart (SP21WC, SP24WC, SP27WC, SP30WC, and SP33WC)



Materials

Shell: **HDPE core with aluminum exterior**

Base: **16-gauge steel with powder-coat finish**

Bumpers: **ABS thermoplastic polymer**

Corner extrusions: **Anodized aluminum**

Work surface: **Antimicrobial ABS thermoplastic polymer**

Casters: **Heavy duty, 5" full-swivel, ball bearing with steel plate and frame**

Dimensions (nominal)

Pace 21

Exterior: **29.25"w x 21.25"d x 33"h**

Weight: **110 lbs**

Pace 24

Exterior: **29.25"w x 21.25"d x 36"h**

Weight: **120 lbs**

Pace 27

Exterior: **29.25"w x 21.25"d x 39.25"h**

Weight: **135 lbs**

Pace 30

Exterior: **29.25"w x 21.25"d x 42.5"h**

Weight: **145 lbs**

Pace 33

Exterior: **29.25"w x 21.25"d x 45.5"h**

Weight: **155 lbs**

Options

- 5 case colors
- Back- and side-mounted accessories
- 3"- and 4"-high casters for under-counter applications
- Any combination of drawer and FlexCell configurations can be specified for your cart as long as the combined height of the drawers and FlexCell sections equals the height of the cart case



Storage solutions for healthcare

1239 Comstock Street | Marne, Michigan 49435
888-435-2256 phone | 616-743-5944 fax | innerspacehealthcare.com